

CL4OVT14

HORIZONTAL FREEZER -40 °C 140 LITERS

[Download PDF](#)

- ✓ temperature range from -45 °C to -10 °C
- ✓ internal volume 140 liters
- ✓ n. 1 cover with self balancing hinges
- ✓ n. 1 lock with key
- ✓ microprocessor electronic control
- ✓ digital display that shows the current temperature (°C/°F)
- ✓ audible and visual alarms
- ✓ n. 2 different alarms: high temperature (adjustable threshold), probes defective
- ✓ evaporator static type
- ✓ manual defrosting
- ✓ default configuration: 1 basket
- ✓ n. 4 turning castors
- ✓ internal structure made in non-toxic plastic material
- ✓ ecological polyurethane insulation 80 mm thick
- ✓ external dimensions LxDxH: 720 x 695 x 890 mm
- ✓ weight 55 kg
- ✓ electrical connection: 1/N/PE 230 V - 50 Hz - 300 W
- ✓ power cord with Schuko plug

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 8211 - fax +39 0522 821453 - service.instruments@smeg.it