

PR377

---



REFRIGERATOR +4 °C 381 LITERS

[Download PDF](#)

---

- ✓ temperature range from +2 °C to +15 °C
- ✓ internal volume 381 liters
- ✓ n. 1 self-closing insulated double glass door
- ✓ n. 1 safety door lock with key
- ✓ microprocessor electronic control
- ✓ digital display that shows the current temperature (°C/°F)
- ✓ audible and visual alarms
- ✓ n. 4 different alarms: door open, temperature too high, temperature too low, probe failure alarm
- ✓ automatic defrosting system
- ✓ evaporator ventilated with forced air
- ✓ shelves made of stainless steel wire protected with non-toxic PVC
- ✓ default configuration: 5 adjustable shelves in height (optional, up to 10)
- ✓ alternative configuration: non toxic plastic drawers with dividers (optional, up to 10)
- ✓ n. 2 back castors + n. 2 front mounting adjustable feet
- ✓ internal structure made in non-toxic plastic material
- ✓ ecological polyurethane insulation 40 mm thick
- ✓ external dimensions LxDxH: 600 x 595 x 1860 mm
- ✓ weight 82 kg
- ✓ electrical connection: 1/N/PE 230 V - 50 Hz - 350 W
- ✓ power cord with Schuko plug

