

DTP140CA



PROFESSIONAL DUAL POSITIVE TEMPERATURE +4 °C 700/700 LITERS

[Download PDF](#)

- ✓ temperature range from +2 °C to +15 °C
- ✓ internal volume 1400 liters (700 + 700)
- ✓ n. 2 blind self-closing doors
- ✓ n. 2 safety door locks with key
- ✓ double microprocessor electronic control
- ✓ double soft touch user interface with LCD display
- ✓ n. 4 languages available (Italian, English, French, Spanish)
- ✓ memory capacity of 4 MB - 46000 data record, last 20 alarms available on display
- ✓ audible and visual alarms
- ✓ n. 7 different alarms: high temperature (adjustable threshold), low temperature (adjustable threshold), open door, condenser cleaning, probes defective, power failure, backup battery failure
- ✓ 12 V backup battery (optional) for temperature visualization and alarms acquisition up to 48h without power connection
- ✓ USB port for data downloading
- ✓ GSM module (optional)
- ✓ Ethernet port for LAN connection (optional)
- ✓ defrosting with automatic evaporation of condensate water (programmable cycles) for both compartments
- ✓ evaporator ventilated with forced air for both compartments
- ✓ illumination made by vertical LED for positive compartment
- ✓ default configuration: n. 8 shelves made in AISI 304 stainless steel and height adjustable both for 1 and 2 compartments (optional, up to 14)
- ✓ alternative configuration: drawers made in AISI 304 stainless steel with Plexiglas frontal side both for 1 and 2 compartments (optional, up to 14)
- ✓ n. 6 wheels + n. 2 front mounting adjustable feet
- ✓ internal structure made in AISI 304 stainless steel
- ✓ ecological polyurethane insulation 60 mm thick
- ✓ external dimensions LxDxH: 1400 x 800 x 1990 mm

- ✓ weight 229 kg
- ✓ electrical connection: 1/N/PE 230 V - 50 Hz - 520 W
- ✓ detachable power cord with Schuko plug

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 8211 - fax +39 0522 821453 - service.instruments@smeg.it

